



Product Code	28042
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279280425

## Description

### CLEOPATRA ELECTRIC CHAIR - 3 motors - light beige

Electrically controlled 3 motors chair to adjust separately height, backrest and seat inclination by foot switch or push button panel.  
Footrest adjustable by gas lift pump.  
Reversible arms and adjustable headrest.  
Cushion included.

Operating power: 100-240 V - 50/60 Hz.

Manual: GB, FR, IT, ES.

Size: 196 x 82 x h 58-100 cm

Net weight: 88 kg

Max. load: 175 kg

Height adjustment: 58-100 cm

Footrest extension: 15 cm adjustable by gas lift pump

Face hole behind cushion