



Product Code	80450
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279804508

## Description

### VET DIGITAL THERMOMETER

- flexible tip
- waterproof
- LCD size: 15 x 7 mm

Multilanguage box and manual: GB, FR, IT, ES, PT, GR.

## Technical Specifications

- Measurement range: 32 ~ 42.9 °C
- Accuracy: ± 0.1 °C
- Dimension: 138x21x10 mm
- Battery: 1,5V LR41
- Measurement time: 60 sec.