

VISIOFOCUS VET NON CONTACT THERMOMETER with Bluetooth



Product Code	80457
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8022386066014

Description

VISIO FOC US® VET - Bluetooth

VisioFocus Vet can measure body temperature of animals at a distance in few second without touching their skin. Thanks to the optical aiming system, temperature can be taken on the eyes through a LED system that indicates the correct distance from the skin. The temperature can also be taken on rectum area, auricle or gum. Built-in cap protects from dust and dirt.

Automatic adjustment to the environments: the software automatically applies a correction factor according to ambient conditions, just leaving it for at least 3-5 minutes in the room. Alternatively you can run a manual calibration (Manual Quick Calibration System, MQCS) in only 3 seconds.

Suitable for all animals

Multipurpose: used to take also the temperature of objects, food, liquids or other surfaces.

Bluetooth system to connect with smartphone: VisioFocus Vet App helps to choose appropriate animal and record a large number of data.

Supplied with 4 AAA batteries.

Technical Specifications

Body measurement range: 32-43.5°C (89.6-110.3°F)
General measurements range: 1-80°C (33.8-176 °F)
Resolution (body - general): 0.1°C (0.1°F)
Dimension: 144x43.5x21.5 mm (5.66x1.71x0.85") including cap
Battery: 4 AAA - LR03
Battery life is approx. 3 years or 30,000 measurements
Measurement time: immediate
Distance from the animal (Calculated by an optical aiming system): about 6 cm (2.4 inches)