



Product Code	28048
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279280487
---------	---------------

Description

SABA ELECTRIC chair - 4 motors - avio blue

Electrically controlled 4 motors chair to adjust separately angle of the backrest, seat inclination, footrest and lifting of the whole chair.

Hand controller and built-in foot controller which allow operator to adjust chair.

High-grade upholstery leather easy to clean and fresh feeling.

Headrest with face hole.

Supplied with removable cushion and remote controller.

Metal base for high stability.

Operating power: 100-240V - 50/60Hz

Manual: GB, FR, IT, ES.

Size: 185 x 88 x h 63-90 cm

Net weight: 98 kg

Max load: 175 kg

Height adjustment of the whole chair 63-90 cm

Angle adjustment of the seat inclination 0° to 20°

Angle adjustment of the backrest 0° to 80°

Footrest angle adjustment 0° to 65°