

RIESTER RI-SCOPE VET OPERATION OTOSCOPE HEAD HL 2.5V - 10541



Product Code	80480
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC 4045396010789

Description

RI-SCOPE® OPERATION VET OTOSCOPE HEAD 2.5V - 10541 needs 80488

- swivel lens with 2.5 X fold magnification
- moveable, slotted speculum holder can be locked in any position
- open design for unhindered use during examination and operation
- bayonet fitting for fast and secure attachment to the handle
- all-metal operation otoscope for long service life
- fast exchange of lamp at the front of the otoscope