

RIESTER VET OTOSCOPE SET HL 2.5V AA-handle - 3853



Product Code	80490
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	4045396153189
---------	---------------

Description

VET HALOGEN OTOSCOPE 2.5V AA handle - 3853

Ri-scope® operation otoscope

Metal specula, unslotted, sterilisable:

- Ø 4 mm, 83 mm long
- Ø 5 mm, 62 mm long
- Ø 6 mm, 82 mm long

Metal specula, slotted, sterilisable:

- Ø 7.5 mm, 68 mm long
- Ø 9.4 mm, 92 mm long

1 spare lamp
swivel lens with 2.5 X magnification