



Product Code 80550

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279805505

Description

Electronic sphygmomanometer for automatic measurement of animals blood pressure. Large, high definition colour LCD display.

- data storage for up to 100 groups of data with date and time
- five cuff sizes: small, small/medium, medium, large, X-large
- automatic one-touch operation
- alarm function with possibility to set alarm limits
- dual scales: kPa and mmHg
- automatic power off
- optional SpO2 probe with clip for tongue or ear

Poor results with small (poor signal) and restless (more artifact) animals.

Supplied with small/medium cuff (10-19 cm) and multilingual manual (GB, FR, IT, ES, PT).

Technical Specifications

Measurement method: oscillometric
 Measurement range: 0-39.6 kPa (0-297 mmHg)
 Resolution: mmHg
 Accuracy: static pressure: ± 0.4 kPa (± 3 mmHg)
 Inflation - Deflation: automatic inflation - automatic multistep deflation
 Display: 2.8" colour LCD screen
 Power: four "AA" batteries or optional power adapter
 Safety: class II device type BF applied part
 Dimension - Weight: 130x110x80 mm - 300 g (without batteries)

SpO2 - optional
 SpO2 range: 0-100% (resolution 1%)
 SpO2 accuracy: 70%-100%: $\pm 2\%$
 PR range: 30bpm-250bpm, (resolution 1bpm)
 PR accuracy: ± 2 bpm or $\pm 2\%$ (select larger)