

# SpO2 PROBE for 80550 - need 80558



Product Code	80557
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	6945040109884

## Description

Optional SpO2 probe for VET electronic sphygmomanometer (code 80550)

SpO2 range: 0-100% (resolution 1%)

SpO2 accuracy: 70%-100%:  $\pm 2\%$

PR range: 30bpm-250bpm, (resolution 1bpm)

PR accuracy:  $\pm 2$ bpm or  $\pm 2\%$  (select larger)