



Product Code	80710
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279807103

Description

GIMA 900 VET MULTIPARAMETER MONITOR

Multi-parameter Vet device which can monitor the vital physiological parameters: non-invasive blood pressure (NIBP) and pulse oxygen saturation (SpO2).

- 3.5"LCD display for remote review
- flexible configurations of SpO2 and NIBP
- history data records reviewed in waveform, list and trend graph
- USB connection for data output and software upgrade

Technical Specifications

NIBP

Measurement mode: Manual, Auto, STAT

Measurement range:

- SYS mmHg 40-265
- MAP mmHg 27-222
- DIA mmHg 20-200

SpO2

Measuring range: 1%~100%

Measuring accuracy: 2% for range 70 to 100%

Pulse Rate

Measuring range: 30~400 bpm

Other Specifications

Power supply: 100~240 VAC, 50/60 Hz

Battery: 11.1 V / 4,400 mAh (Li-ion Battery)

Size: 168x172xh 242 mm

Standard accessories

- VET SpO2 probe
- NIBP VET Suntech cuff set (Set of 5 cuffs: 3-6 cm, 4-8 cm, 6-11 cm, 7-13 cm, 8-15 cm)
- Lithium rechargeable battery
- Power supply cable
- Dust protective mantle
- Manual (GB, IT)