



Product Code	80740
Unit of sale	1 kit
Minimum order	1
Type	No medical device
Class	

UK-REP	Yes
CH-REP	Yes
EAN/UPC	8023279807400

Description

CRYOMEGA® VET DUAL DELIVERY CRYOSURGERY KIT

CryOmega® Vet Dual Delivery is the perfect choice for cryosurgery treatment of benign lesions.

Cryogen reaches a temperature of -68 C° and is immediately ready to use.

It can be delivered by cones or buds.

It combines the portability of a canister with the effectiveness found in more expensive liquid nitrogen-based units.

Single handed design, easy to use, no sutures required.

Multilingual manual: GB, FR, IT, ES, DE.

Available GR, RO, NL, BG, SE, PL, PT from www.gimaitaly.com

- Effective freezing
- Anesthesia not required
- High efficacy
- Dual delivery: buds or cones
- Cones: Different sizes to provide precise treatment
- Buds: Appropriate when treating hard to reach or thin-skin surfaces

Technical Specifications

Type of lesion treatments / Freeze time range

Warts: 3 x 10 seconds

Eyelid tumors: 10 - 20 seconds

Meibomian gland tumors: 10 - 20 seconds

Skin tags: 2 x 5 seconds

Distichia: 10 seconds

Papillomas: 3 x 10 seconds

Adenomas: 3 x 10 seconds

Lick granulomas: 3 x 10 seconds

Anal growths: 3 x 10 seconds

CRYOSURGERY ADVANTAGES

Less time involved

Tissue sparing

Less bleeding

Less post-surgical pain

Less sedation or none required

Minimal cutting or none required

Standard accessories

1 canister of cryogen gas
handle and trigger
4 cone adapters
70 cones (3mm, 5mm, 7mm, 9mm, 12mm)
1 cone selection guide
20 buds
1 neoprene practice pad