



Product Code	80805
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279808056

Description

SP-20 VET PULSE OXIMETER

SP-20 handheld pulse oximeter can simultaneously measure SpO2 and Pulse Rate by placing sensor on tongue, ear, leg or tail. Colour TFT display shows plethysmogram, Perfusion Index and bar-graph.

- Lightweight, small and easy to carry
- Continuous monitoring
- Up to 500 hours data storage for SpO2 and PR
- Data management software for analysis
- Audible and visual alert function
- Vertical and horizontal display

SP-20 pulse oximeter is supplied with universal "Y" type SpO2 probe and 2 adapters: big and small clip for different measuring, orange rubber cover, charging base and user manual (GB, IT).

Technical Specifications

SpO2 range: 1100% \pm 2%
 Pulse rate range: 30400 bpm \pm 2%
 Perfusion Index range: 0.2~20%
 Display: 3.5" colour TFT display
 Size: 158x73x25 mm
 Net weight: 230 g
 Power supply: 2000 mAh lithium battery