



Product Code	28097
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	237553
CND	Z12159002
EAN/UPC	8023279280975
PARAF	927604470
GTIN	08023279280975
GMDN	35457

## Description

### EOLO NEBULIZER - piston

Innovative, portable, easy to use nebulizer for home use, recommended to nebulize antibiotics and bronchodilator drugs. Plastic body with high thermal insulation. Supplied with adult and pediatric face mask, bulb, mouthpiece, air tubing, air filter. Multilingual user manual and box (GB, FR, IT, DE, ES, PT, GR, PL\*, Arabic\*). Warranty 2 years.

\*Box only, not manual

Made in Italy.

## Technical Specifications

Operating voltage: 230V - 50/60Hz - 170VA / 110V or 60Hz on request  
 Operating pressure: 1.1 bar  
 Air flow: 14 l/min  
 Nebulization rate 0.4 ml/min  
 Particle size < 5 µm  
 Max Pressure: 2.5 bar  
 Dimensions: 22x14.5x12 cm - 1.7 kg  
 Noise level (EN 13544-1) 55 db  
 Class of risk IIa  
 Norms: insulation class. Type B