



Product Code	28102
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1348898
CND	Z12159002
EAN/UPC	8023279281026
PARAF	909847687
GTIN	08023279281026
GMDN	35457

Description

MISTRAL NEBULIZER - piston

Compact portable nebulizer for home use. Duty cycle: 20 min on / 40 min off. Supplied with Hi Flo bulb, adult and pediatric mask with string, mouthpiece, nosepiece, air tubing and air filter. Multilingual user manual and box: GB, FR, IT, ES, PT, DE.

Made in Italy.

Technical Specifications

Operating voltage: 230V - 50/60Hz - 170VA
Operating pressure: 1.1 bar
Air flow: 14 l/min
Nebulization rate 0.35 ml/min
Particle size < 5 µm
Max Pressure: 2.5 bar
Dimensions: 13x8.9x14.6 cm - 1.5 kg
Noise level (EN 13544-1) 60 db
Class of risk IIa
Norms: insulation class. Type B