


JAR HOLDER 2 I for 28194-5,28118,28200,28203 needs 28109 - spare



| | |
|--------------|-------|
| Product Code | 28108 |
|--------------|-------|

| | |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

| | |
|---------------|---|
| Minimum order | 1 |
|---------------|---|

| | |
|------|------------|
| Type | Spare part |
|------|------------|

| | |
|-------|--|
| Class | |
|-------|--|

| | |
|--------|----|
| UK-REP | No |
|--------|----|

| | |
|--------|----|
| CH-REP | No |
|--------|----|

| | |
|---------|---------------|
| EAN/UPC | 8054619108468 |
|---------|---------------|

Description

Jar holder 2 I for 28194-5, 28118, 28200, 28203

Needs 28109