



Product Code	28183
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2059697
CND	Z120105
EAN/UPC	8023279281835
PARAF	990144192
GTIN	08023279281835
GMDN	36777

Description

Supervega evo battery aspirator

Supervega 118 and Supervega Evo are electric medical and portable devices for nasal, oral and tracheal aspiration of body fluids in children or adults.

Two or three options of operation: by 100-240 V power feeding, 12 V car adapter or internal rechargeable battery.

Equipped with a long-life brushless motor eliminating any type of smell and coal residuals.

Large LCD

Clear reading of vacuum values along with soft keys for vacuum adjustment increase the accuracy of aspiration.

FEEDBACK system

It controls and manages the power of aspiration, providing a long autonomy to the lightweight lithium-ion battery and quite noise level during operation.

PROXIMITY function: NO-TOUCH ON/OFF

This function prevents and avoids cross-contamination inbetween patients.

Various options of jar

Standard autoclavable (121°C) Makrolon® polycarbonate, optional FLOVAC® disposable liners or autoclavable (143°C) APEC® polycarbonate.

In compliance with EN 60601-1-11 for use in the home healthcare environment and carries IP22 degree protection.

Made in Italy.

Technical Specifications

- Motor: Oil-less and maintenance-free piston pump
- Power:
 - Internal rechargeable Lithium battery (14.8 V - 5.2 A) - Autonomy 70 min
 - Car-ambulance adapter (12 V - 4 A)
 - AC/DC adapter 110-240 V, 50-60 Hz
- Max vacuum (adjustable): -0.75 Bar; -75 kPa; -563 mmHg
- Max free air flow rate: 26 l/m
- Size: 35x19x15 cm Weight: 2.7 kg
- Duty cycle: non stop operation (by mains or 12 V only)
- Noise level min/max: 46.7/61.8 dBA

Standard accessories

Polycarbonate jar 1,000 ml with overflow valve system

Antibacterial filter

Suction catheter

Silicone tube set with conical connector

Universal power adaptor (100-240V - 50/60 Hz)

12V car adapter

User Manual GB, FR, IT, DE, ES