

"SUPER VEGA ON TROLLEY" SUCTION UNIT



Product Code	28191
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	590729
CND	Z120105
EAN/UPC	8023279281910
PARAF	987472875
GTIN	08023279281910
GMDN	36777

Description

SUPERVEGA on trolley - suction 2x2 l jar

Continuously adjustable vacuum, with indicator.
Silicone connection tube 8x14 mm. 2 jars 2 l, made of sterilizable polycarbonate (autoclavable at 121°C).
Approved for NON-STOP operation without over heating.
Suction power: 40 air litres/min.
Steel painted frame, fitted with 5 star base.

Vacuum indicator Max suction: 600 mm Hg - 0.80 bar
Multilingual user manual: GB, FR, IT, DE, ES.

Supervega on trolley also fits Flovac® jars 28272 for disposable liners 28273.

Made in Italy.

Technical Specifications

Size: 32x30xh 90 cm Weight 6.2 kg
Operating voltage: 220-230 V - 50/60 Hz; other on request
Power consumption: 110 VA
Maximum suction: -0.80 bar (600 mm Hg)
Flow: 40 air litres/minute
Operating cycle: continuous with AC/DC
Noise level: 61.5 dBA
Norms CEI 62-5 (IEC 601-1) - 93/42 EEC - Class II - Type B