

"SUPER VEGA ON TROLLEY" SUCTION UNIT with footswitch



Product Code	28192
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	590729
CND	Z120105
EAN/UPC	8023279281927
PARAF	987472887
GTIN	08023279281927
GMDN	36777

Description

SUPERVEGA on trolley - suction 2x2 l jar with footswitch

Continuously adjustable vacuum, with indicator.

Silicone connection tube 8x14 mm. 2 jars 2 l, made of sterilizable polycarbonate (autoclavable at 121°C).

Approved for NON-STOP operation without over heating.

Suction power: 40 air litres/min.

Steel painted frame, fitted with 5 star base.

Vacuum indicator Max suction: 600 mm Hg - 0.80 bar

Multilingual user manual: GB, FR, IT, DE, ES.

Supervega on trolley also fits Flovac® jars 28272 for disposable liners 28273.

Made in Italy.

Technical Specifications

Size: 32x30x90 cm Weight 6.2 kg

Operating voltage: 220-230 V - 50/60 Hz; other on request

Power consumption: 110 VA

Maximum suction: -0.80 bar (600 mm Hg)

Flow: 40 air litres/minute

Operating cycle: continuous with AC/DC

Noise level: 61.5 dBA

Norms CEI 62-5 (IEC 601-1) - 93/42 EEC - Class II - Type B