

# "SUPER VEGA BATTERY ON TROLLEY" SUCTION UNIT



Product Code	28193
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2059700
CND	Z120105
EAN/UPC	8023279281934
PARAF	987472901
GTIN	08023279281934
GMDN	36777

## Description

Supervega Battery on trolley - suction 2x2 l jar

Professional portable suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children.  
Suction power of 36 l/ min.

Double collection jars on a stand with five castors for easy transport.

Two different working options: AC/DC adapter or rechargeable battery make it a versatile suction unit approved for NON-STOP operation without overheating.

Visual and acoustic alarm indicating low battery level, vacuum-meter (bar & kPa) and vacuum regulator.

User manual: GB, FR, DE, IT, ES.

Made in Italy.

## Technical Specifications

Operating voltage: 100-240 V - 50/60 Hz  
Internal rechargeable Pb battery: 12V - 4A  
Battery autonomy: 60 min  
Recharging: 240 min  
Operating cycle: continuous with AC/DC  
Flow: 36 air litres/minute  
Noise level: 65.5 dBA  
Size - Weight: 32x30xh 99 cm - 7.9 kg