

"MANUALE TOBI" SUCTION ASPIRATOR



Product Code	28220
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	590729
CND	Z120105
EAN/UPC	8023279282207
PARAF	903149829
GTIN	08023279282207
GMDN	36616

Description

TOBI MANUALE - manual suction aspirator

Independent of any power source, Tobi Manuale is simple to operate pressing pedal by foot, hand or knee.

It is used to clear patient's airways safely and efficiently, anytime, anywhere. Precision and solidly built to assure long service life.

Polycarbonate suction bottle with capacity 400 cc, autoclavable at 120°C and spill-over protection.

Bacterial filter included.

Multilingual user manual: GB, FR, IT, DE, ES, PT.

Made in Italy.

Technical Specifications

Weight: 1.15 kg
Size: 225x165x85 mm
Max suction -0.40 bar
Suction power: 22 l/min
Bottle capacity: 0.4 l