


BOTTLE 4 l WITH COVER - autoclavable 121°C



Product Code	28228
Unit of sale	1 pc
Minimum order	1
Type	Spare part
Class	

EAN/UPC	8023279282283
---------	---------------

Description

Jar 4 l with cover - Makrolon®
- Autoclavable up to 121°C
- Compatible with 4 l aspirators.

Autoclavable jars with overflow valve system
Equipped with screw cover with handle for easy grip, autoclavable silicone o-ring for tight seal, overflow valve system integrated in the cover, clear graduated scale in ml with 100 or 200 ml intervals.
Suitable for central vacuum systems and Gima suction units.
Available in two types of polycarbonate, Makrolon® and Apec®.