



Product Code	28243
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	590729
CND	Z120105
EAN/UPC	8023279282436
PARAF	987472863
GTIN	08023279282436
GMDN	36777

Description

SUPER VEGA BATTERY - plastic case - 16 l/min.

Professional, portable suction units for body, oral, nasal or tracheal liquids aspiration in adults or children. Supplied with three different operation options: AC/DC adapter, internal lead rechargeable battery and 12 V make it a versatile unit approved for NON-STOP operation without overheating. Visual and acoustic alarm indicating low battery level. Ergonomic handle ensures a comfortable and secure grip. Bottle with safety float control valve to prevent overflow.

Made in Italy

Technical Specifications

Suction power: 16 air l/min
 Max suction: -0.75 bar (563 mm/Hg)
 Operating Voltage: 100-240 V - 50/60 Hz
 Battery autonomy: approx 80 min
 Recharging time: 240 minutes
 Noise level: 63 dBA
 Size - Weight: 350x210x180 mm - 3.5 kg
 Norms: Class II Type BF. EN 60601-1-1. EN 60601-1-2. EN 60601-1-11. 93/42 EEC

Standard accessories

Bottle 1 l - autoclavable at 121°C
 AC/DC universal adapter (100-240V - 50-60Hz)
 Silicone tubing set
 Anti-bacterial filter
 Suction catheter
 Cigarette lighter adaptor (12V 4A)
 User manual (GB, IT, FR, DE, ES, GR, PL, RO, PT)