

ASPEED SUCTION ASPIRATOR - 230V single pump



Product Code	28244
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1758504
CND	Z120105
EAN/UPC	8033717010750
PARAF	987471253
GMDN	36777

Description

Piston-type continuous cycle electric aspirators give high performance and great durability.
Equipped with a protective thermal cut-out relay.
They require no maintenance or lubrication.
A motor-protector cap totally prevents aspirated liquids or secretions from reaching and damaging the vacuum pump.
Ideal for clinical use.
Made in Italy.

- Aspeed 15 l
- Power: 230 V
- Pump: single
- Case: metal

Technical Specifications

Operating voltage: 230 V-50 Hz
Bottle capacity: 1 l
High vacuum: high flow
Air flow: 15 l/min
Adjustable vacuum level: 0÷ -0.85 bar (0÷ -85 kPa)
Weight: 3.5 kg
Case material: metal
Noise level: 53 dBA

Standard accessories

Autoclavable polycarbonate jar 1,000 cc with safety valve (overflow protection)
Disposable suction liner 1 l
Antibacterial hydrophobic filter
Sterile disposable cannula
Sterile manual flow regulator
Set of atoxic sterilizable silicone tubes
Power Cable
User Manual: GB, FR, IT, DE, ES