



Product Code	28300
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	364316
CND	Z120606
EAN/UPC	8019781472038
UMDNS	12415
PARAF	921453243
GMDN	35169

Description

MAG 2000 PEMF MAGNETOTHERAPY - 2 channels - 300 Gauss

Ideal devices for domiciliary and professional magnetotherapy. 35 programmes: 20 pre-set programmes for the treatment of several pathologies, 14 programmes identified by the working frequency (from 1 to 100 Hertz) and an Autoscan programme with continuous frequency variation for the treatment of both hard and soft tissues at the same time. MAG 2000 PLUS has also the possibility to store 10 user memories, choosing the frequency up to 120 Hz, intensity and therapy time: a great possibility for the professional user who need to use a specific treatment for each patient. Possibility to set the intensity of magnetic field for all the programmes up to 300 Gauss (150 Gauss each channel) for MAG 2000 or up to 400 Gauss (200 Gauss each channel) for MAG 2000 Plus to easily treat different body parts through elastic therapeutic belt with 3 solenoids. Supplied with elastic therapeutic belt with 3 solenoids, power supply, carrying bag and manual (GB, FR, IT, ES, DE).

Size: 180x110x50 mm

Technical Specifications

35 PROGRAMMES

20 pre-set programmes:

- Osteoporosis
- Arthrosis
- Arthritis
- Cervical arthrosis
- Articular pain
- Sprains
- Fractures
- Epicondylitis
- Epitrocleitis
- Intercostal Contusion
- Lumbalgy
- Lumbar pain
- Shoulder arthrosis
- Knee arthrosis
- Periarthritis
- Coxarthrosis
- Muscular Atrophy
- Muscular contracture
- Osteonecrosis

14 working freq. + autoscan:

- 1 Hz
- 3 Hz
- 5 Hz
- 10 Hz
- 15 Hz
- 20 Hz
- 30 Hz
- 40 Hz
- 50 Hz
- 60 Hz
- 70 Hz
- 80 Hz
- 90 Hz
- 100 Hz
- autoscan