

IONO BASE PLUS - 2 channels



Product Code	28306
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No
RDM (NSIS)	177135
CND	Z120601
EAN/UPC	8023279283068
PARAF	987472255
GMDN	35372

Description

IONO BASE - PLUS - 2 channels

Local administration of medicines by Ionophoretics, electric osmosis and treatment of hyper hydrosis (sweating disease) and cellulitis. Current generator for hydro galvanization, aesthetic shaping and against cellulitis slight blemish. Microprocessor controlled, with selection of parameters and internal programs. Simple and safe to use, although it maintains characteristics of a professional equipment.

Supplied with ABS case, velcro elastic bands, rubber conductive electrodes 80x120 mm and 50x50 mm and envelopes sponges for electrodes.

Made in Italy.

Technical Specifications

Operating Voltage: 230 V - 50/60 Hz
Size: 138x68x193 mm
Absorbed Max Power: 10 VA
Max Output Current: 20 mA
Safety class: II BF - 93/42/EEC.