



|               |                |
|---------------|----------------|
| Product Code  | 28309          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2586237       |
| CND        | Z120610       |
| EAN/UPC    | 8019781472175 |
| UMDNS      | 11248         |
| PARAF      | 921453294     |

Description

MIO-SONIC ULTRASOUND THERAPY

Innovative and high performance device both for home use and professional use.

It combines compactness and high performances in a small , light and easy to use applicator. The lightness and handiness allows a continuous use without any difficulty.

Mio Sonic can be used for the following applications:

- pain therapy
- sport injuries
- aesthetic treatments

Supplied with gel bottle, carrying and storage bag, medical power supply and multilanguage manual: GB, FR, IT, ES. On request DE.

Technical Specifications

9 programs;  
1 MHz ultrasound, maximum output power 6.4 W, maximum intensity 1.6 W/cm2;  
Collimated ultrasound beam for an optimal penetration beam;  
3 ultrasound beam intensity levels: 0,08 W/cm2 (low) - 0,8 W/cm2 (medium) -1,6 W/cm2 (high);  
Power supply: 100-240V, 50/60 Hz  
Direct contact device (not suitable for water immersion)  
Size-weight: 200x47x69 mm - 190 g