

MAG 3000 MAGNETOTHERAPY - 2 channels



Product Code	28314
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2369981
CND	Z120606
EAN/UPC	8019781401083
PARAF	985798103

Description

MAG 3000 PEMF MAGNETOTHERAPY - 2 channels

MAG 3000 is a certified medical device for pulsed electromagnetic field magnetotherapy (PEMF). Suggested for the treatment of delayed union, osteoporosis, bone oedema, osteonecrosis, as well as ulcers and neuropathies. Thanks to the high intensity of the magnetic field, it is also suitable for the treatment of bone fractures even in the presence of rigid bandages or plaster casts.

- two independent channels
- colour display
- magnetic field intensity up to 150 gauss per channel
- frequency from 1 to 100 Hz
- 41 programmes available
- possibility to manage two treatments at the same time by setting different times and programmes (double user)
- list of the last 10 programmes used

For professional users (doctors, therapist,...) and for in-home therapy.

Operating power 110-240 V.

Technical Specifications

41 PROGRAMMES:

- Osteoporosis
- Knee arthrosis
- Healing
- Arthrosis
- Scapulo humeral periarthritis
- Cutaneous ulcers
- Arthritis
- Coxarthrosis
- Discopathy
- Cervical arthrosis
- Muscular atrophy
- Myalgia
- Articular pain
- Muscular treatment
- Neuropathy
- Cervicalgias
- Osteonecrosis
- Muscle strain
- Chronic pain
- Algodystrophy

Muscular cramp
Fractures
Cartilage lesion
Rizarthrosis
Epicondylitis
Ligament lesion
Impingement syndrome
Epitrocleitis
Bone cedema
Carpal tunnel
Pseudoarthrosis
Tendinitis
Titanium prosthesis
Lumbalgy
Chondropathy
Rotator cuff
Lumbar pain
Antinflammatory
Tarsal tunnel
Shoulder arthrosis
Whiplash syndrome

Standard accessories

3 solenoids belt applicator
Medical power supply
User manual (GB, FR, IT, ES) on request DE.
Magnet to verify therapy operation
Carrying bag
Three non-woven fabric strips 15 x 150 cm