

# ELASTIC BELT 40x8 cm for reusable electrodes



|               |                |
|---------------|----------------|
| Product Code  | 28330          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2432639        |
| CND        | C0280          |
| EAN/UPC    | 8033003264966  |
| GTIN       | 08033003264966 |

## Description

### UNIVERSAL ELASTIC BELTS FOR REUSABLE ELECTRODES

Elastic belts - blue - 40x8 cm

Elastic belts for reusable electrodes are easy to use and assure a very good seal.  
Made of synthetic fibres (nylon, elastomer and transparent nylon yarn) to avoid skin irritation in patients.  
Belts are available in 4 different lengths to match almost all kinds of reusable electrodes on the market.