



|               |                |
|---------------|----------------|
| Product Code  | 28338          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2586120       |
| CND        | Z120606       |
| EAN/UPC    | 8019781472090 |
| PARAF      | 921487132     |
| GMDN       | 47593         |

Description

MAG 1000 MAGNETOTHERAPY - 2 channels - 20 programmes

Mag 1000 is a high frequency PEMF device particularly recommended for the treatment of pathologies affecting soft tissues like muscles and tendons. It offers 20 pre-set programmes with 3 adjustable intensity levels. The device generates an electric signal with a carrier frequency of 27 MHz, which is modulated with specific frequencies (from 8 Hz to 640 Hz) and therefore is able to soothe the anatomical site affected by pain. It is projected and indicated for healing, rehabilitation and functional recovery treatments of many pathologies. Continuous operating mode. It is intended for both professional and homecare use.

Size: 153x91x38 mm

Supplied with elastic therapeutic belt, carrying bag, rechargeable battery pack 2,300mAh, battery charger, user manual in GB, FR, IT, ES, DE.

Technical Specifications

- 20 programmes:
- Arthritis-Arthrosis
  - Bruises
  - Rheumatisms
  - Myalgia
  - Pain-articular prosthesis
  - Tendonitis
  - Herniated disc
  - Hematoma
  - Muscular pain
  - Epicondylitis
  - Osteoporosis
  - Epytrocleitis
  - Fractures
  - Lumbago
  - Sprains
  - Muscular contractures
  - Dislocations 2h
  - Medium application
  - Traumas 4h
  - Long application