

TAP2000 MATTRESS FOR MAGNETOTHERAPY



Product Code	28339
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2307002
CND	Z12069099
EAN/UPC	8019781101525
UMDNS	12415
PARAF	975062579
GMDN	35169

Description

TAP 2000 MAGNETOTHERAPY MAT for MAG 2000

TAP 2000 is a 40×40 cm mat with 2 cables of 3 solenoids each (total 6 solenoids) ideal for magnetotherapy treatments. It is compatible with high intensity/low frequency pulsed electromagnetic field devices: MAG 2000, MAG 2000 Plus and MAG 2000 Premium.

Supplied with one connection cable and user manual.