



Product Code	28370
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2581709
CND	Z120601
EAN/UPC	8019781471109
UMDNS	13762
PARAF	921453128
GMDN	35372

Description

MIO-CARE TENS - 20 medical programmes

20 PROGRAMMES
300 APPLICATIONS

- 2 independent channels stimulator
- compensated biphasic square wave: avoids hazardous thermic effect of polarization
- easy to use, digital keyboard
- remaining therapy time and programmes instruction shown on graphic display
- internal rechargeable battery pack
- maximum intensity 100 mA (200 mA pk-pk)
- compliance with 93/42/EEC MDD

Supplied in a carrying bag with 2 connection cables with splitters, user manual* (GB, FR, IT, DE, ES), 8 adhesive pads, internal rechargeable battery with battery charger. *PT on request

Dimensions: 60x135x19 mm

Technical Specifications

20 MEDICAL PROGRAMMES:

- Tens conventional (rapid)
- Tens endorphinic (delayed)
- Tens maximum values
- Inflammatory
- Neck pain/cervicogenic headache
- Backache/sciatic pain
- Sprains/bruises
- Vascularization
- Muscle relaxant
- Hand and wrist pain
- Plantar stimulation
- Epicondylitis
- Epitroclea
- Periarthritis
- Haematomas
- Atrophy prevention
- Atrophy (trophism rehabilitation)

