


GELLED ELECTRODES 45x35 mm



Product Code	28387
Unit of sale	bag of 4 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2269232
CND	N010299
EAN/UPC	8033695723789
PARAF	977812294
GMDN	17191

Description

DISPOSABLE GELLED ELECTRODES - latex free
High material flexibility allows easy application of electrodes on treated areas. Solid gel maintains the electrode adhesion and conductivity unchanged during procedure.