



Product Code	28405
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2986637
CND	Z12062899
EAN/UPC	8023279284058
PARAF	986873457

Description

GIMA-CARE TENS - 60 medical programmes

Dual channel multi-functional electrotherapeutic device generates electric stimulation to help easing aches and pain, muscle training and relaxation.

Easy to use, Gima-care has 2 pre-set treatment programmes with adjustable treatment time and 1 custom programme for each body part, allowing the user to adjust treatment time and pulse rate and width.

- 3 in 1 Tens + Ems + Massage
- professional tens waveform (variable pulse frequency and width)
- body part illustration
- 60 programmes (10 body parts, 40 levels)
- powerful output (max. 120 mA) (at 500 ohm load)
- solid and clear design

Continuous operation for around 10 treatments if fully charged.

Supplied with 4 pads 50x50 mm, 2 output cables, 1 USB cable.

User manual (GB, FR, IT, ES, DE).

Technical Specifications

60 programmes - TENS/EMS/MESSAGE:

TENS (Transcutaneous Electrical Nerve Stimulation)
Body parts: neck, shoulder, arm, hand, back, abd, hip, leg, foot, joint
Treating programme: 30 programmes, 10 adjustable programmes

EMS (Electrical Muscle Stimulation)
Body parts: neck, shoulder, arm, hand, back, abd, hip, leg, foot
Treating programme: 27 programmes, 9 adjustable programmes

MESSAGE
Methods: knead, rub, tap
Treating programme: 3 massage manipulations to all body parts

TECHNICAL SPECIFICATIONS:
Channel: 2
Waveform: Bi-phase square-wave pulse
Output current: Max. 120 mA (at 500 ohm load)
Output intensity: 0 to 40 levels, adjustable

Power supply: 3.7 V Li-ion rechargeable batteries

Display: backlight white display

Dimension - Weight: 109x54.4x23 mm - 85g

Classification: BF type applied part, internal power equipment, IP22

Norms: IEC60601, IEC61000-4-2, IEC61000-4-4, IEC61000-4-5, IEC61000-4-11, IEC61000-4-8, IEC61000-4-3