


GEL AIR 2D ANTIDECUBITUS CUSHION 41x41x7.5 cm



Product Code	28538
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2044549
CND	Y033303
PARAF	975385903
GTIN	03701127718379

Description

GEL AIR 2D CUSHION 41x41x7.5 cm

PREVENTION FOR HIGH TO VERY HIGH DECUBITUS RISK

Load 30-120 Kg.

Made in France.

Moulded viscoelastic foam base (i.e. Visco Mouss) with a polyurethane gel insert. The viscoelastic fluid gel (blue pad) ensures maximum pressure relief of vulnerable bony areas while the regular gel and the foam increase the contact surface and maximize pressure relief. The cushion's shape ensures a good posture and prevents from sliding forward. Foam density: 86 kg/m³

GIMA line of seat cushions is designed to prevent decubitus risk for patients on wheelchair or permanently on chairs. All cushions in Gel or Foam are supplied with a waterproof PU cover with zipper and handle.