

COMPRESSOR for 28564



Product Code	28565
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2398246
CND	V080701
EAN/UPC	8023279285659
PARAF	975386006

Description

Special silent pump for an excellent pressure regulation by a knob on the front panel.

It can be hung to any kind of hospital bed frame by means of two adjustable hooks.

For use with bubble or alternating cell mattresses.

Technical Specifications

Size: 25x12xh 10 cm
Operating voltage: 110-220 V, 50-60 Hz
Working pressure: 40-100 mmHg
Air output: 6-7 liter/min
Cycle time: 6 min
Case: TUV-CE IEC/EN60601-1; 60601-1-2