

INTERCHANGEABLE CELL AIR MATTRESS - stage II



Product Code	28567
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2398254
CND	V080701
EAN/UPC	8023279285673
PARAF	975386018

Description

INTERCHANGEABLE CELL-MATTRESS

Recommended for the prevention and treatment of pressure sore stage I and II in medium-term and domiciliary therapies.

The 20 9.5 cm-high Nylon+PVC cells provide superior comfort.

Cells (20+1 spare) are fixed by means of pressure stud and can easily be replaced in case of breakage.

Manual and box in 10 languages: GB, FR, IT, ES, DE, PT, PL, HU, GR, Arabic.

Technical Specifications

Total size: 200 x 86 x 9.5 cm
Cell height: 9.5 cm (4")
Material: Nylon+PVC
Thickness: 0.32 mm
Suitable for persons up to: 140 kg