

COMPRESSOR for 28567



Product Code	28568
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2398246
CND	V080701
EAN/UPC	8023279285680
PARAF	975386020

Description

- Special silent pump for an excellent pressure regulation by a knob on the front panel.
- It can be hung to any kind of hospital bed frame by means of two adjustable hooks.
- For use with bubble or alternating cell mattresses.

Technical Specifications

Size: 25x12xh 10 cm
Operating voltage: 110-220 V, 50-60 Hz
Working pressure: 40-100 mmHg
Air output: 6-7 liter/min
Cycle time: 12 min
Case: TUV-CE IEC/EN60601-1; 60601-1-2