

# KIT 8050 ALTERNATING PRESSURE MATTRESS + PUMP



Product Code	28578
Unit of sale	1 kit
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	817830
CND	Y033399
EAN/UPC	6944749500541
PARAF	975386071

## Description

Advanced systems of alternating pressure mattress.

The multi-mode pump offers.

- alternating mode: the pressure is periodically redistributed by inflating and deflating cells in order to mimic natural sleep movement patterns and avoid long term pressurization of tissue.
- constant static mode: the pressure is about 2/3 lowered when compared to the same pressure level from alternating mode by means of immersion and envelopment. Visual alarms notify the caregiver unusual situation, reducing the risk of human error.

Suitable for the prevention and treatment of all stages of pressure ulcers.  
Manual in: GB, IT.

## Technical Specifications

### Mattress:

- Total size: 200 x 90 x 22 cm
- Alternating: 1 in 2 (group A and B)
- Cell height: 22 cm (8,66")
- Cover material: Nylon PU with cotton lining
- Thickness: 0.3 mm
- Cells material: TPU
- Suitable for persons up to: 180 kg
- Number of cells: 20 + 1 (spare)

### Pump:

- Size: 30 x 13.5 x h 11 cm
- Cell height: 12.5 cm
- Working pressure: 55-120 mmHg
- Operating voltage: 220 V, 50 Hz
- Air output: 7-8 liter/min
- Cycle time: 12 min
- Case material: Flame retardant ABS - TUV-CE IEC/EN60601-1; 60601-1-2
- Static Mode, low pressure alarm, quick connector, transfer function