

INFRARED LAMP 250 W - desk



Product Code	28652
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8436562860097
PARAF	977812395

Description

INFRARED LAMP - 250 W - desk with table clamp

Infrared lamp with heat distributed evenly. The luminaire is mounted with a highly precise articulated arm. The arm allows easy adjustment in different positions, providing a long reach and maximum vertical and horizontal movement. Metallic structure. Lamp is manufactured with durable materials and with very low maintenance costs.

Head Ø 21 cm
Arm length 90 cm (45+45)
Internal spring balanced arm allows different positions

Made in Spain.