



Product Code	28653
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8436562860097
---------	---------------

## Description

### INFRARED LAMP - 250 W - trolley

Infrared lamp with heat distributed evenly. The luminaire is mounted with a highly precise articulated arm. The arm allows easy adjustment in different positions, providing a long reach and maximum vertical and horizontal movement. Metallic structure. Lamp is manufactured with durable materials and with very low maintenance costs.

Head Ø 21 cm

Arm length 90 cm (45+45)

Internal spring balanced arm allows different positions

Made in Spain.