

WAIST WRAP HEATING PAD



Product Code	28675
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279286755
PARAF	950454153

Description

WAIST WRAP HEATING PAD

Waist heating pad made of super cozy fleece surface with stamp welding manufacturing. 3 temperature settings with level indicator. Automatic switch off after 90 minutes. Detachable hand controller allows to wash the pad. Fast heating and overheating protection.

Supplied with multilingual box and manual: GB, FR, IT, DE, NL, PT, GR, SE, PL, RO, FI, CZ.

Technical Specifications

Size: 29x69 cm
Power consumption: 100 W
Voltage: 220-240 V