

TRULENE NON ABSORB. SUTURE gauge 2/0 circle 3/8
needle 24 mm - 70 cm - blue



Product Code	22872
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781302
CND	H010201010403
EAN/UPC	8023279228724
GMDN	13909

Description

TRULENE - monofilament non absorbable suture - polypropylene
Equivalent Ethicon brand: PROLENE

- Gauge: 2/0
- Metric: 3.0
- Colour: Blue
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 24 mm
- Suture length: 70 cm

Technical details:

- monofilament polypropylene
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- hernia repair
- permanent fixtures
- general surgery
- micro surgery
- orthopedic surgery
- cardiovascular surgery
- plastic surgery
- can also be used in ophthalmic and neurological procedures

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.