


KIT MEN'S EXTERNAL CATHETER diam 35 mm + HYDROCOLLOID STRIP



Product Code	28725
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2491339
CND	U070201
EAN/UPC	7290018953082
PARAF	950080337

Description

KIT: MALE EXTERNAL CATHETER + ADHESIVE STRIP IN A POUCH

External male catheters made from natural medical latex.

The graded cap design eliminates the possibility of twisting of the drainage system and blockage of the urine flow.

No leaking or ballooning is possible.

- Kit: catheter Ø 35 mm + hydrocolloid strip
- Size: XL

Multilingual polybag and instructions (GB, FR, IT, ES, PT, DE, NL, PL, CZ, RO, HU, SE, FI, GR, BG)