

BD PLASTIPAK SYRINGES WITHOUT NEEDLE - 30 ml Luer Lock



Product Code	23873
Unit of sale	box of 60 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1869585
CND	A020102020102
GTIN	30382903012290
GMDN	47017

Description

BD PLASTIPAK SYRINGES WHITOUT NEEDLES - LUER LOCK

- Capacity: 30 ml
- Model: Luer - lock
- BD Code: 301229