

# PEHA 991078 STANDARD ANATOMIC FORCEPS - straight - 14 cm



Product Code	38749
Unit of sale	box of 40 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Not required
CH-REP	Not required
RDM (NSIS)	2795862
CND	V9012
EAN/UPC	4052199761596

## Description

### STERILE FORCEPS

- Standard anatomic forceps straight - 14 cm - 991078 - box of 40

Single-use instruments made of matt brushed, high-quality stainless steel.

Made in EU, they offer similar characteristics as multi-use instruments in terms of intended performance, haptics and design.

Sterile, individually packed and ready-to-use for a variety of surgical interventions.

They allow precise handling during procedures. Peha® color lock® safety marking allows quick and reliable differentiation from multi-use instruments.

The self-adhesive traceability label on the packaging facilitates easy, fast and less error prone data entry into the patient records.

Excellent choice for health care facilities with no processing unit.

Recyclable and disposable in existing clinical waste containers after use.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, GR, HU, CZ, SI, RU, Arabic.