

TRULENE NON ABSORB. SUTURE gauge 4/0 circle 3/8
needle 19 mm - 75 cm - blue



Product Code	22879
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	1781304
CND	H010201010403
EAN/UPC	8023279228793
GMDN	13909

Description

TRULENE - monofilament non absorbable suture - polypropylene
Equivalent Ethicon brand: PROLENE

- Gauge: 4/0
- Metric: 1.5
- Colour: Blue
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 15 mm
- Suture length: 75 cm

Technical details:

- monofilament polypropylene
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- hernia repair
- permanent fixtures
- general surgery
- micro surgery
- orthopedic surgery
- cardiovascular surgery
- plastic surgery
- can also be used in ophthalmic and neurological procedures

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.