

SMEDLEY HAND DYNAMOMETER



Product Code	28791
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279287912
PARAF	977812698

Description

SMEDLEY HAND DYNAMOMETER - aluminium

Spring type hand dynamometer used to measure grip strength.
Adjustable handle comfortably accommodates any hand size.
Dial indicator remains at maximum until reset and is calibrated both in kg (0-100) and pounds (0-220).
Manual in: GB, FR, IT, ES.

Approximate measurement of values.