

HYPOR FILM DRESSING FOR FIXATION 10x25cm



Product Code	58843
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2514624
CND	M04010202
EAN/UPC	8015143008808
PARAF	987475365

Description

HYPOR FILM DRESSING FOR FIXATION - STERILE

- Hypor film 10 x 25 cm

Sterile adhesive dressing on transparent polyurethane film apts to fix catheters, gauzes, bandages and other dressings.
The films acts as water and bacterial barrier.
Permeable to air and water vapor, lets the skin breathe and avoids the risk of maceration. Individually packed.
Extremely elastic substrate and highly conformable to sensitive areas.
Reduces the risk of fluid leaking and infection.

Multilingual box: GB, FR, IT, ES, DE.
Medical Device Class I.

Rounded edges to prevent the dressing from peeling off.
Hypoallergenic acrylic adhesive, dermatologically tested.
Protective liners easy removable with gloves.
Made in Italy.