

STERILE SPLINTER PLAIN FORCEPS - 12.5 cm



Product Code	38885
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	1847836
CND	K01020104
EAN/UPC	8023279388855
GTIN	08023279388855
GMDN	62479

Description

- Splinter plain forceps - 12.5 cm

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.

Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.