

MULTIPLE KEY TURNER



Product Code	28941
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279289411
PARAF	977812799

Description

Designed to assist person with limited wrist strength or movement.
People with arthritis or other limited hand mobility can manipulate keys with ease.
House keys and most car keys (without plastic head covers) fit into the slots and are held by a single screw.

Can hold 1 to 3 keys. Keys fold into handle when not in use to slip into pocket or purse.

Length 13 cm.