

HANDLE GRIP CUP



Product Code	28967
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279289671
PARAF	951853391

Description

HANDLE GRIP CUP

Transparent plastic resistant mug with two large handles to allow a more secure grip.
Extrawide base for stability and to reduce tipping over.
It includes a spouted lid with a little hole for straws.
Suitable for patients with limited hand control and tremors.
Mug only is dishwasher and microwave safe.

Capacity: 295 ml
Diameter: 7.60 cm (3")
Height: 12.7 cm (5")
Weight: 0.27 kg